



Raya Hyaluron Gel 2.5

INCI breakdown Aqua, Alcohol, Sodium Hyaluronate, Glycerin, Sodium Levulinate, Sodium Anisate, Citric Acid

Trade name: **Raya Hyaluron Gel 2.5**

Application: Cosmetics

Packaging: in cans, bottles and drums

Filling quantity: 1 kg/unit
5 kg/unit
Special packaging upon request

Preservatives: Ethanol,
p-anisic acid/levulinic acid

Shelf life: up to 36 months,
not opened

Storage conditions: cool (10 – 20 °C)
protected from light and heat



Product Specification Hyaluronic acid Gel 2.5

Test	Specification	Result
Product code:		CEX 52 M02
Appearance	clear liquid	conforms
Dynamic Viscosity 25 °C Haake Viskotester Air	700 – 900 mPa.s	780 mPa.s
Density (20 °C) (Ph.Eur.2.2.5)	1.013 – 1.017 g/cm ³	1.015 g/cm ³
Assay Hyaluronic Acid	2.3 – 2.7 %	2.5 %

Used Hyaluronic Acids:

Raya Hyaluron	SHMW	(1200 – 1800 k Da)
Raya Hyaluron	MMW	(200 – 400 k Da)
Raya Hyaluron	LMW	(10 – 50 k Da)
Raya Hyaluron	Oligomer	(< 10 k Da)